

Prosthetic Care Instructions

Care of the Silicone Liner

The silicone liner should be cleaned every evening with one of the following techniques.

Once per week:

- Turn the liner inside out.
- Spray with rubbing alcohol.
- Wipe clean with a lint free towel.
- Turn liner right side out and hang or lay flat to dry.

Daily:

- Turn the liner inside out.
- Wash with warm water and a small amount (dime size) of mild soap.
- Do not use Dial antibacterial soap!**
- Rinse thoroughly.
- Wipe clean with a lint free towel.
- Turn liner right side out and hang or lay flat to dry.

Care of Prosthetic Socks

Prosthetic socks should be washed each day after use.
Most socks can be machine washed in warm water and tumbled dry.
Nylons can be machine washed and hung to dry.
Follow the instructions printed on the package from the manufacturer.

Care of the Prosthesis and Foam Liner

The socket and foam liner of the prosthesis should be sprayed with rubbing alcohol and wiped clean weekly.

Care of the Prosthetic Foot

The sock and shoe should be removed and inspected once weekly to check for any debris that may be causing excessive wear on the prosthetic foot.

If you have any questions regarding the care or fit of the prosthesis, please contact your prosthetist.