

Cosmetic Restoration

A cosmetic restoration is done by professional artists. This is the most cosmetic (attractive) and most expensive covering available. A mold is taken of the amputee's sound side limb; then a custom made cover is fabricated and painted to match the sound side limb.

For more information visit Alternative Prosthetic Services at www.alternativeprosthetics.com.

Indications:

- Can be applied at almost any level of involvement and sometimes in conjunction with a functional orthosis/prosthesis.



Photos courtesy of Prosthetic Laboratories



Modern Technology. Traditional Values. Proven Techniques.