

## **Hip Knee Ankle Foot Orthosis (HKAFO)**

The molded plastic hip knee ankle foot orthosis is a custom fabricated thermoplastic orthosis with metal uprights. It extends from the hip to the foot to control motion and alignment of the knee and ankle joints as well as mediolateral (side to side) stability of the hips. There are different components used to lock the knee if needed. Joints may also be added to the ankle if motion is desired.

### **Indications:**

- Spina bifida, muscular dystrophy, traumatic paraplegia, cerebral palsy, osteogenesis imperfecta, and rotational issues at the hip.



*Modern Technology. Traditional Values. Proven Techniques.*