

Leather Ankle Foot Orthosis (AFO)

The leather ankle foot orthosis is a custom orthosis with padded plastic reinforcement that is covered with leather on the inside to protect the skin, as well as on the outside for durability. This orthosis provides mild mediolateral (side to side) ankle stability.

Indications:

- Arthritis, severe pes planus with pain, posterior tibial tendon dysfunction, Charcot arthropathy, and degeneration of the foot/ankle.



Photo courtesy of Prosthetic Laboratories



Modern Technology. Traditional Values. Proven Techniques.