

Metal and Leather Knee Ankle Foot Orthosis (KAFO)

The metal and leather knee ankle foot orthosis is custom fabricated and extends from the thigh to the foot. It is designed to control motion and alignment of the knee and ankle. As with the plastic KAFO, different components can be used to lock the knee and allow motion at the ankle if desired. It is designed for minimal skin coverage and is attached to a shoe.

Indications:

- Weak quadriceps, knee contractures, genu varum/valgum deformities, post-operative stabilization, genu recurvatum, femoral fractures, post ligament repair, spina bifida, cerebral palsy, muscular dystrophy, paraplegia, and post-polio.



Photo courtesy of Prosthetic Laboratories



Modern Technology. Traditional Values. Proven Techniques.