

Plastic Knee Ankle Foot Orthosis (KAFO)

The molded plastic knee ankle foot orthosis is a custom fabricated thermoplastic orthosis with metal uprights. It extends from the thigh to the foot to control motion and alignment of the knee and ankle. There are different components used to lock the knee if needed. Joints may also be added to the ankle if motion is desired.

Indications:

- Weak quadriceps, knee contractures, genu varum/valgum deformities, post-operative stabilization, genu recurvatum, femoral fractures, post ligament repair, spina bifida, cerebral palsy, muscular dystrophy, paraplegia, and post-polio.



Photo courtesy of Cascade