

Denis-Browne Splint

The Denis-Browne splint consists of a spreader bar with foot plates that attach to the bottom of the patient's shoes. It is used to keep the patient's feet in the desired amount of internal/external rotation and control curvature of the foot (club foot).

Indications:

- Talipes equinovarus (clubfoot) and calcaneovalgus, congenital convex pes planovalgus, metatarsus adductus, and external femoral/tibial pathologic conditions.



Photo courtesy of Prosthetic Laboratories

PROSTHETIC
LABORATORIES

Modern Technology. Traditional Values. Proven Techniques.