

## Pediatric Hip Abduction Orthosis

The pediatric hip abduction orthosis is partially flexible and foam-lined. The brace maintains desired hip abduction, but is designed to allow walking. It is commonly used following the “Pavlik Harness” as a full-time or night-time only intervention.

### Indications:

- Unilateral or bilateral developmental hip dysplasia in children aged six months to three years.



Photo courtesy of Prosthetic Laboratories

**PROSTHETIC**  
**LABORATORIES**

*Modern Technology. Traditional Values. Proven Techniques.*