

## **TLSO for Scoliosis**

The TLSO for Scoliosis is a custom made plastic orthosis that is worn to control curvatures of the spine. This orthosis is primarily used to correct curves from T8 and lower. A mold of the back and pelvis are done on a special frame while corrective forces are applied. The plastic orthosis is fabricated to exact dimensions using x-rays, measurements and the mold of the spine.

### **Indications:**

- Idiopathic/congenital scoliosis, kyphosis, lordosis, and neuromuscular deformities at or below the level of T8.



Photo courtesy of Prosthetic Laboratories



*Modern Technology. Traditional Values. Proven Techniques.*