

## **BK Custom Shaped Cover and Flexible Protective Outer Surface Covering System**

### **Custom Shaped Cover**

A custom shaped cover is flexible foam which is attached to the outside of a definitive prosthesis and then shaped to resemble the sound side limb as closely as possible.

### **Flexible Protective Outer Surface Covering System**

The flexible protective outer surface covering protects the prosthesis from excessive wear and adds cosmesis (makes it more attractive). This “skin” is applied to the foam once the patient receives a definitive prosthesis and the custom shaped cover is complete.



Photo courtesy of Prosthetic Laboratories



*Modern Technology. Traditional Values. Proven Techniques.*