

AK Custom Shaped Cover and Flexible Protective Outer Surface Covering System

Custom Shaped Cover

A custom shaped cover is made of flexible foam that is attached to the outside of a definitive prosthesis and then shaped to resemble the sound side limb as closely as possible.

Flexible Protective Outer Surface Covering System

A flexible protective outer surface covering protects the prosthesis from excessive wear and adds cosmesis (makes it more attractive). This “skin” is applied to the foam once the patient receives a definitive prosthesis and the custom shaped cover is complete.



Photo courtesy of
Prosthetic Laboratories



Modern Technology. Traditional Values. Proven Techniques.